

Labexchange

The World's No.1



Your TopAir Partner in Germany!

www.labexchange.com

Your partner

Discover our services as a certified partner for all Top Air products. With a comprehensive range of top-quality equipment always in stock, we ensure swift access to the solutions you need.

Get in touch:

Hubert Sauter

✉ hubert.sauter@labexchange.com

☎ +49 (0)7475-9514-24

Biosafety Cabinet



WORLD PREMIERE:
Integrated Particles
Monitoring System!

TopAir's Class II A2 Biological Safety Cabinet protects lab staff, the environment and sensitive work processes in which biological agents are applied.

Features:

- Color touch screen.
- Two ULPA H15 filters @99.9995% @ 0.1 um.
- Maintenance & technical faults alarms.
- Timers and counters management .
- Germicidal water proof UV light system.
- Airflow Pattern: 70% circulation, 30% ex-haus.

Ductless Fume Hood



Topair's Polypropylene Ductless Fume Hood provides a safe work environment for lab staff working with acids and harsh chemicals.

Features:

- Air velocity display.
- Color touch screen.
- 10-speed fan.
- Hour counter for filter.
- Filter replacement alarm for HEPA/carbon filters.
- Variable Air Volume System.
- Chemical Sensor Alarm.

Fume Cupboard



The Polypropylene Fume Hoods protect laboratory staff from noxious fumes released by acids, dangerous gas and organic solutions.

Features:

- Excellent chemical resistance.
- Color touch screen.
- Air suction from both the top and back panel.
- Air velocity monitoring and display.
- Optional: sink/water tap/gas tap/vacuum tap, electric sockets.



TOPAIR
CLEAN AIR SOLUTIONS

Labexchange – Die Laborgerätebörse GmbH
Bruckstr. 58, 72393 Burladingen, Germany

www.labexchange.com